

Georgia Adopt-A-Stream Bacterial Monitoring




*2 Martin Luther King Jr. Drive
Suite 1462, East Tower
Atlanta, Georgia 30334
AdoptAStream.Georgia.gov
404.463.1464*

Georgia Adopt-A-Stream

- What is it?
 - Georgia's volunteer water quality monitoring program
- Program Goals
 - **A**: Increase public **awareness**
 - **D**: Collect quality baseline water quality **data**
 - **O**: Gather **observations**
 - **P**: Encourage **partnerships** between citizens & local government
 - **T**: Provide **tools & training**

Bacterial Monitoring

- Involves: collecting, plating, and counting E. coli colonies
- Purpose: **quickly assess health risks due to bacterial contamination of surface waters**
- Characterizes: Stream health by identifying potential for harmful bacteria in stream



EPA Quality Assurance Project Plan

- Quality Assurance/
Quality Control (QA/QC)
- Only individuals are certified
- Certification is valid for one year
- Volunteers must attend an annual recertification workshop
- Only certified volunteers can submit data



To Become a Certified QA/QC Volunteer

FIELD & LAB:

Volunteers must demonstrate how to properly collect and plate a sample

WRITTEN TEST:

Volunteers must pass a written evaluation with a score of at least 80%

&

Volunteers must correctly identify E.coli colonies and calculate E.coli levels of example plates with accuracy of at least 90%

What is a Watershed?

- A watershed is the land area from which water, sediment, and dissolved materials drain to a common point along a stream, wetland, lake or river.
- Its boundaries are marked by the highest points of land around the waterbody.



When, Where, & How Often?

Where: Same site location and in a well mixed area of flowing water

When: Same time of day and during normal flow conditions

How often: At least **once a month**



Safety Considerations


If conditions are too dangerous to sample...

DON'T SAMPLE!

- Wait until storm has stopped and strong flow has subsided
- Never sample alone
- Remember to wear gloves and boots
- Disinfect and dispose of used plates properly
- Receive permission from land owner before going onto private property

What are Bacteria?


- **Single-celled, living microscopic organisms**
- There are more bacteria on Earth than any other living thing
- Found all over the world in many different environments



Why are Bacteria Important?

Benefits

- Decomposition
 - Digestion
- Nutrient cycling
- Can be used to aid in pollution control
 - Sewage treatment
 - Oil spills
 - Other pollutants




Risks

- Some bacteria that live in soil, water or air can cause human, animal and plant health problems
 - Pathogenic bacteria


Coliform Bacteria

- Total Coliform refers to a biological family of bacteria that are naturally found in soil
- Fecal coliforms are a subgroup within coliform bacteria found in intestinal tracts of humans and other warm-blooded animals
- E. coli are one subgroup of fecal coliform



Why Monitor for E. coli ?

- E. coli serve as an indicator species
 - **High levels indicate possible presence of pathogens**
- Sources of E. coli in our waterways
 - **Wildlife**
 - **Livestock**
 - **Urban storm runoff**
 - **Leaking pipes**
 - **Failing septic systems**



Weather and Seasonal Influences




- **Bacteria numbers often increase following a heavy rainstorm or heavy runoff event**
- E. coli counts are often higher in summer compared to lower counts in winter months

Warmer surface water temps = higher replication rates for E. coli

Colder surface water temps = lower replication rates for E. coli

6 Steps for E.coli Monitoring

1. Preparing the blank/control
2. Collecting a sample
3. Plating your samples and blank
4. Incubating
5. Counting
6. Disposing



1. Preparing the blank/control

Creating a blank checks to see if sampling methods allow for contamination.

- Label Whirl-pak® bag as a blank
- Wear gloves and remove the perforated seal
- Use small white tabs to pull open the bag
- Fill the bag 2/3 full with distilled water (brought from home)
- Whirl!
- Place blank in a sanitized cooler with ice where other sample will be stored




2. *Collecting a sample*

- Use one Whirl-pak[®] to sample upstream of yourself
- Clearly label your sample
- **Place sample in your cooler after collection**
 - **Exposure to UV will decrease bacteria levels**
- **Plate sample as soon as possible, but the maximum holding time on ice is no more than 24 hours**



3. *Plating your sample*

- Clean work area with disinfectant spray
- Before plating, invert sample to mix the water
- Prepare 3M Petrifilm plates
 - Check expiration dates!
 - Label: Blank, 1, 2, 3 (+ site name)
- **Run in triplicate + 1 blank (3 plates from your one sample + 1 control/blank plate)**
 - **4 plates total**



4. Incubating

- **35°C ± 1 degree for 24 hours ± 1 hour**
 - Turn on early to get stable temperature reading with a digital thermometer inside
 - Adjust using metal dial on top
- Check minimum and maximum temperatures after incubating




5. Reading the Results...

counting the *E.coli* colonies


- Standard reporting units for bacteria are in **CFU/100 ml → Colony Forming Units per 100ml**
- Only count **blue colonies with entrapped gas bubbles!**
- Do not count colonies that are growing more than halfway off of the medium

Possible gas bubble patterns associated with gas producing colonies.

All of these examples would be counted when reading plates.



The Blank/Control




There should not be any colonies on the blank.

If any colonies appear on blank, sample is null and void! And new sample must be taken from site location.


Example #1

How many *E. coli* colonies can you see?




Example #2

How many *E. coli* colonies can you see?



Example #3


How many *E. coli* colonies can you see?



Example #4


How many *E. coli* colonies can you see?

Do not count blue colonies without gas bubbles




Example #5

How many *E. coli* colonies can you see?




Example #6

How many *E. coli* colonies can you see?




Example #7

How many *E. coli* colonies can you see?




Example #8

How many *E. coli* colonies can you see?



Example #9




How many *E. coli* colonies can you see?

Would this plate be acceptable as a blank?

Example #10

How many *E. coli* colonies can you see?




Example #11

How many *E. coli* colonies can you see?


When plates turn a deep red color, or plates have many gas bubbles due to the presence of many coliforms, count ALL blue colonies as presumptive *E. coli* colonies

This is the only case where blue colonies without gas bubbles are counted



Example #12

How many *E. coli* colonies can you see?



How Do We Calculate Results from the Plates?




Plate 1




Plate 2




Plate 3

	Plate 1	Plate 2	Plate 3
<i>E. coli</i> <i>Colonies</i>	3	2	3

How Do We Calculate Results?

Calculate the results:

$$\text{Step 1: } \frac{\text{Plate 1} + \text{Plate 2} + \text{Plate 3}}{3 \text{ plates}} = \text{Average CFU/1 ml}$$

$$\text{Step 2: } \text{Average CFU/1 mL} \times 100 \text{ mL} = \# \text{ CFU/100 ml}$$

How Do We Calculate Results?

	Plate 1	Plate 2	Plate 3
<i>E. coli</i> Colonies	3	2	3

Step 1: $\frac{3 + 2 + 3}{3} = 2.67 \text{ CFU/1 ml}$

Step 2: $2.67 \text{ CFU/1 ml} \times 100 \text{ ml} = \underline{\underline{267 \text{ CFU/100 ml}}}$

How Do We Calculate Results from these Plates?




Plate 1




Plate 2




Plate 3

	Plate 1	Plate 2	Plate 3
<i>E. coli</i> Colonies	TNTC	TNTC	TNTC

How Do We Calculate Results?

	Plate 1	Plate 2	Plate 3
<i>E. coli</i> Colonies	TNTC	TNTC	TNTC

Step 1:
$$\frac{\text{TNTC} + \text{TNTC} + \text{TNTC}}{3} = \text{TNTC}$$

*TNTC is generally defined as more than 150 colonies

6. Disposal and Clean-up

- **Spray plates with disinfectant (10% bleach, disinfectant spray, rubbing alcohol, etc.) and seal in bag or used Whirl-Pak before throwing away**
- Spray incubator with disinfectant to clean after use
- Remember to wear gloves – these are living bacteria colonies!
- Wash your hands when finished

Storage of Petrifilm

- If using within one month, keep in the fridge
- If not, keep them in the freezer for storage
- Remember to thaw before use!



EPA Recommended E. coli Levels for recreational waters




	<i>Designated swimming</i>	<i>Moderate swimming area</i>	<i>Light swimming area</i>	<i>Infrequent swimming area</i>
<i>E. coli (cfu/100ml)</i>	<i><235</i>	<i><298</i>	<i><410</i>	<i><576</i>

These levels correspond to an acceptable risk level of 8 people out of 1000 getting sick
(from US EPA 1986, 2002a)


Fecal vs. E. coli

- According to the EPA, approximately 60% of fecal coliform represents *E. coli*
 - For example:
 - A fecal coliform sample that has 1000 cfu/100 ml has approximately 600 cfu/100 ml of *E. coli*
 - A fecal coliform sample that has 500 cfu/100 ml has approximately 300 cfu/100 ml of *E. coli*

But what level is AAS concerned with?



Counts that Warrant Action



GEORGIA ADOPT-A-STREAM: Bacterial Form

To be conducted every month

All monitoring programs

Bacterial specific

SITE INFORMATION	Group Name: _____		Event Date: _____ (MMDDYYYY)	
	Group ID: G-_____ Site ID: S-_____		Time Sample Collected: _____ (HHMM am/pm)	
WEATHER	Stream Name: _____		Time Spent Sampling: _____ (Min)	
	Monitor(s): _____		Total Time Spent Traveling (optional): _____ (Min)	
	Number of Participants: _____		Furthest Distance Traveled (optional): _____ (Miles)	
OBSERVATIONS	Present conditions (check all that apply)		Amount of rain, if known?	
	<input type="checkbox"/> Heavy Rain <input type="checkbox"/> Steady Rain <input type="checkbox"/> Intermittent Rain <input type="checkbox"/> Overcast <input type="checkbox"/> Partly Cloudy <input type="checkbox"/> Clear/Sunny		Amount in Inches: _____ In Last Hours/Days: _____ *Refer to wunderground.com for rainfall data	
	Flow/Water Level: (check all that apply) <input type="checkbox"/> Dry <input type="checkbox"/> Stagnant/Still <input type="checkbox"/> Low <input type="checkbox"/> Normal <input type="checkbox"/> High <input type="checkbox"/> Flow (over banks)			
	Water Clarity: <input type="checkbox"/> Clear/Transparent <input type="checkbox"/> Cloudy/Somewhat Turbid <input type="checkbox"/> Opaque/Turbid			
	Water Color: <input type="checkbox"/> No Color <input type="checkbox"/> Brown/Muddy <input type="checkbox"/> Green <input type="checkbox"/> Milky/White <input type="checkbox"/> Tannic <input type="checkbox"/> Other: _____			
	Water Surface: <input type="checkbox"/> Clear <input type="checkbox"/> Oily Sheen: does it break when disturbed? Yes/No (circle one) <input type="checkbox"/> Algae <input type="checkbox"/> Foam <input type="radio"/> Greater than 3" high <input type="radio"/> It is white			
	Water Odor: <input type="checkbox"/> Natural/None <input type="checkbox"/> Gasoline <input type="checkbox"/> Sewage <input type="checkbox"/> Rotten Egg <input type="checkbox"/> Fishy <input type="checkbox"/> Chlorine <input type="checkbox"/> Other: _____			
	Photos: Please take images to document your observations and changes in water quality conditions. Photo point directions can be found in the manuals. Images can be submitted online with your other data.			
	Trash: <input type="checkbox"/> None <input type="checkbox"/> Yes, I did a cleanup <input type="checkbox"/> This site needs an organized cleanup			
	BACTERIAL	3M Petrifilm Method: <i>Escherichia coli</i>		
Run three (3) plates/tests for each site, plus one (1) blank plate. Process within 6-24hrs, incubate at 35°C ±1° and read at 24 ± 1 hr				
Plate		Colonies	Find AVG of Number of Colonies	cfu/100mL
Blank			(total # colonies/total # of plates (do not include blank)	
1			(/) x 100 =	
2		Sample Holding Time (HH): _____		
3		Date START (MMDDYYYY): _____	Date END (MMDDYYYY): _____	
Total # Colonies		Time START (HHMM): _____	Time END (HHMM): _____	
		MIN Temp (°C): _____	MAX Temp (°C): _____	
COMMENTS	Any changes since you last sampled at this site? If yes, please describe.			

Observations

- Flow/Water Level
- Water Clarity
- Water Color
- Water Surface
- Water Odor
- Photos
- Trash



OBSERVATIONS	Flow/Water Level: <input type="checkbox"/> Dry <input type="checkbox"/> Stagnant/Still <input type="checkbox"/> Low <input type="checkbox"/> Normal <input type="checkbox"/> High <input type="checkbox"/> Flood (over banks) <small>(check all that apply)</small>
	Water Clarity: <input type="checkbox"/> Clear/Transparent <input type="checkbox"/> Cloudy/Somewhat Turbid <input type="checkbox"/> Opaque/Turbid
	Water Color: <input type="checkbox"/> No Color <input type="checkbox"/> Brown/Muddy <input type="checkbox"/> Green <input type="checkbox"/> Milky/White <input type="checkbox"/> Tannic <input type="checkbox"/> Other: _____
	Water Surface: <input type="checkbox"/> Clear <input type="checkbox"/> Oily sheen: Does it break when disturbed? Yes/No (circle one) <input type="checkbox"/> Algae <input type="checkbox"/> Foam <input type="radio"/> Greater than 3" high <input type="radio"/> It is pure white <input type="checkbox"/> Other: _____
	Water Odor: <input type="checkbox"/> Natural/None <input type="checkbox"/> Gasoline <input type="checkbox"/> Sewage <input type="checkbox"/> Rotten Egg <input type="checkbox"/> Fishy <input type="checkbox"/> Chlorine <input type="checkbox"/> Other: _____
	Photos: Please take images to document your observations and changes in water quality conditions. Photo point directions can be found in the manuals. Images can be submitted online with your other data.
	Trash: <input type="checkbox"/> None <input type="checkbox"/> Yes, I did a cleanup <input type="checkbox"/> This site needs an organized cleanup

Bacterial Data Form

- Use Bacterial data form (Chemical/Bacterial combo data form may also be used)
- Fill in number of colonies for each plate
 - Blank count should be zero
- Volunteers can send scanned color images of plates to AAS or their local coordinator if they need assistance with counting colonies

BACTERIAL	3M Petrifilm Method: <i>Escherichia coli</i>				
	Run three (3) plates/tests for each site, plus one (1) blank plate. Process within 6-24hrs, incubate at 35°C ±1° and read at 24 ± 1 hr				
	Plate	Colonies	Find AVG of Number of Colonies		cfu/100mL
	Blank		(total # colonies/total # of plates (do not include blank)		
	1		(/) x 100 =		
	2		Sample Holding Time (HH): _____		
	3		Date START(MMDDYYYY): _____	Date END (MMDDYYYY): _____	
	Total # Colonies		Time START (HHMM): _____	Time END (HHMM): _____	
		MIN Temp (°C): _____	MAX Temp (°C): _____		

Submit the Data

As soon as possible after monitoring is complete

Data should be submitted to the state program's **online database:**

AdoptAStream.Georgia.gov

Share your data with partners, local governments and your local Adopt-A-Stream coordinators



Georgia Adopt-A-Stream

Georgia's Volunteer Water Quality Monitoring Program

Search this site



Get Involved

Confluence



Citizen Monitoring

Data Views

Data Entry

Materials & Resources

View: All active sites



Data Submission Form

Register Group or Site

Trainers: Enter Workshop Data

Trainers: Certificates & Letters

Trainer Workshop History

Events shown in time zone: Eastern Time



Check out our most recent newsletter!



Announcements

AAS Volunteer Monitoring Conference - Confluence 2017

Learn More About Getting Started With Adopt-A-Stream

From the website's Home Page, select "Data Submission Form" under the Data Entry tab.

[Site](#) | [Chemical](#) | [Bacterial](#) | [Macroinvertebrate](#) | [Stream Habitat Survey](#)

GEORGIA ADOPT-A-STREAM Data Submission Form

[Trainings calendar](#)[Errors and Warnings list](#)

You must enter Site information and click "submit" at the bottom of the page before moving on to the chemical, bacterial, macroinvertebrate, or stream habitat survey forms. You must click submit on each page on which you enter data.

Below six parameters required

AAS monitors, Total participants, Site, Event Date, Event Time, Time Spent Sampling

You cannot submit a form that has **Errors** or missing **Required Data**.

You can submit a form that has **Warnings**, but it will be flagged as out of compliance with the AAS quality assurance plan.

Site, Weather, and Observations

Site Information

Site:

Enter the site name or site number without the S-, and select from the list.
Note that you must be a member of a group before you can submit data for its sites.

Event date: mm/dd/yyyy <input type="text"/>	Time sample collected: 03 :02 PM hh:mm am/pm	Total number of participants: <input type="text"/>	Time spent sampling: <input type="text"/> minutes	Total time spent traveling: <input type="text"/> minutes <i>Optional</i>	Furthest distance traveled: <input type="text"/> miles <i>Optional</i>
---------------------------------------------------	----------------------------------------------------	-------------------------------------------------------	------------------------------------------------------	--------------------------------------------------------------------------------	------------------------------------------------------------------------------

Participants

Adopt-A-Stream monitors

Enter one at a time, and select from the drop-down list.

Other participants

Enter your site information as well as any weather and observation information on this page.



Georgia Adopt-A-Stream

[Return to Home Page](#)

Georgia's Volunteer Water Quality Monitoring Program

User: **Jennings**

[Citizen Monitoring](#) | [Data Views](#) | [Data Entry](#) | [Reports](#) | [Outreach Staff](#) | [My Profile](#)

[Site](#) | [Chemical](#) | [Bacterial](#) | [Macroinvertebrate](#) | [Stream Habitat Survey](#)

GEORGIA ADOPT-A-STREAM Data Submission Form

[Trainings calendar](#)

[Errors and Warnings list](#)

You must enter Site information and click "submit" at the bottom of the page before moving on to the chemical, bacterial, macroinvertebrate, or stream habitat survey forms. You must click submit on each page on which you enter data.

Below six parameters required

AAS monitors, Total participants, Site, Event Date, Event Time, Time Spent Sampling

You cannot submit a form that has **Errors** or missing **Required Data**.

You can submit a form that has **Warnings**, but it will be flagged as out of compliance with the AAS quality assurance plan.

Site, Weather, and Observations

Site Information

Site:

(scroll to the bottom of the page...)

Comments

Any changes to note since you last sampled at this site?
If so, please describe. Otherwise, please leave blank.

Large text area for entering comments.

Email

Clear check box if you don't want email confirmation.

Click "Submit" at the bottom of the page to record your data. You must submit your site data before you can enter bacterial data



[Site](#) [Chemical](#) [Bacterial](#) [Macroinvertebrate](#) [Stream Habitat Survey](#)

GEORGIA ADOPT-A-STREAM Data Submission Form

[Trainings calendar](#)

[Errors and Warnings list](#)

You must enter Site information and click "submit" at the bottom of the page before moving on to the chemical, bacterial, macroinvertebrate, or stream habitat survey forms. You must click submit on each page on which you enter data.

Below six parameters required

AAS monitors, Total participants, Site, Event Date, Event Time, Time Spent Sampling

You cannot submit a form that has **Errors** or missing **Required Data**.

You can submit a form that has **Warnings**, but it will be flagged as out of compliance with the AAS quality assurance plan.

Site, Weather, and Observations

Site Information

Site:

Enter the site name or site number without the S-, and select from the list.
Note that you must be a member of a group before you can submit data for its sites.

Event date: mm/dd/yyyy <input type="text"/>	Time sample collected: 03 :02 PM hh:mm am/pm	Total number of participants: <input type="text"/>	Time spent sampling: <input type="text"/> minutes	Total time spent traveling: <input type="text"/> minutes Optional	Furthest distance traveled: <input type="text"/> miles Optional
---------------------------------------------------	----------------------------------------------------	-------------------------------------------------------	------------------------------------------------------	-------------------------------------------------------------------------	-----------------------------------------------------------------------


Participants

Adopt-A-Stream monitors

Enter one at a time, and select from the drop-down list.

Other participants

After clicking "Submit," click on the Bacterial tab to continue entering data



[Site](#) | [Chemical](#) | [Bacterial](#) | [Macroinvertebrate](#) | [Stream Habitat Survey](#)

Bacterial Data

3M Petrifilm Method: *Escherichia coli*

Run three (3) plates/tests for each site, plus one (1) blank plate. Process within 0-24 hrs, incubate at 35°C ± 1°, and read at 24 ± 1 hr

Plate					Find AVG of number of colonies total # colonies / total # of plates (do not include blank)	cfu /100 mL
blank	1	2	3	4		
<input type="text"/>	<input type="text"/>	<input type="text"/>	<input type="text"/>	<input type="text"/>	<input type="text"/>	

Sample Holding Time	START		END		TEMPERATURE		Total Time
	Date	Time	Date	Time	Minimum	Maximum	
	<input type="text"/> mm/dd/yyyy	--Select-- hh:mm am/pm	<input type="text"/> mm/dd/yyyy	--Select-- hh:mm am/pm	<input type="text"/> °C	<input type="text"/> °C	

Did you use a method other than 3M Petrifilm Plates? Yes No


Comments

*Any bacterial changes to note since you last sampled at this site?
If so, please describe. Otherwise, please leave blank.*

Fill out the form and click "Submit" to record your bacterial data!

Note: You must click "Submit" on each tab if you are entering data for multiple parameters

Volunteer Monitoring Data Uses



Just the Facts

- *Raise Awareness*
- *Water quality Data*
- *Gather Observations*
- *Encourage Partnerships*
- *Provide Tools and Training*

Quickly assess health risks due to bacterial contamination of surface waters

Single-celled, living microorganisms

CFU/100 mL – colony forming units/100 mL

1000 cfu/100 ml

- *Wildlife*
- *Livestock*
- *Urban storm runoff*
- *Leaking pipes*
- *Failing septic systems*

1 blank + triplicate of sample = 4 plates

Purpose of a blank?

Once a month

Place samples in cooler with ice

Rainstorm

UV light will reduce bacteria levels

Spray plate with disinfectant and seal in plastic bag for disposal

Warmer temps = higher E. coli survival rates

Cooler temps = lower E. coli survival rates

***General coliforms** are red in color with entrapped gas bubbles, while **E. coli** are blue in color with entrapped gas bubbles*

asap, but no more than 24 hours after collection

35°C ± 1 degree for 24 hours ± 1 hour